



Linkfar[®] Pressure Activated Safety Lancet

SAFE:

- All needles are EO-sterilized and protected with a sterility tab.
- Fully enclosed needle housing ensuring that needle is hidden before and after use to prevent sharps injuries.
- Self-destructing mechanism preventing reuse of the device.

COMFORTABLE:

- Silicone-coated, ultra-sharp needles, ideally positioned during skin penetration.
- Various, color-coded product versions to meet specific blood sample applications and address the variety of patient skin types.
- Wide range of designs providing simple, comfortable and secure puncture.

SIMPLE:











- Easy to handle, intuitive activation in 2 steps only with no pre-loading.
- Wide and long safety cap for easy removal – simply twist and pull to remove the needle cover.
- Precise designs for convenient finger positioning and ideal blood sample collection.

Instructions for use

1. Twist off the protective cap.
2. Push the needle into fingertip to get blood sample.
3. Push the cap into the safety lancet back and dispose of the used one to a trash bin.



Product Selection Chart

Item	Color	Specification	Application	Package
04-2122		High flow 21G Needle/2.2mm 	Glucose,Cholesterol,HGB, PT/INR and other like tests	100pcs/box 10,000pcs/case
04-2322		High flow 23G Needle/2.2mm 	Blood glucose testing, Cholesterol,HbA1c.	
04-2618		Medium flow 26G Needle/1.8mm 	Blood glucose testing, HbA1c.	
04-2818		Low flow 28G Needle/1.8mm 	Blood glucose testing, HbA1c.	
04-3015		Low flow 30G Needle/1.5mm 	Pediatric and Frequent testing.	

 **Linkfar**™

Shandong Lianfa Medical Plastic Products Co.,Ltd

Add: No1,Shuangshan Sanjian Road 250200 Zhangqiu City, Jinan, Shandong,China.

Web: www.sdliafa.com

Tel/Fax: 86-531-61328777

Mobile: +86-15098983602

Email: henry@lianyingchina.com